



Pronto Selections
Ready to serve gourmet

Date: Nov 26, 2020

Product Specification

Name: Mac & Cheese

Brand: Pronto Selections

Item #: 8084

UPC: 628451380849

GTIN: 10628451380846



Product Description: Individual macaroni & cheese in a dual ovenable bowl

Product Characteristics

Colour: Yellow/orange

Odour: Cheese

Shape: N/A

Texture: Soft noodles, creamy sauce

Taste: Cheddar

Size/Count: 298g

Product Packaging:

Count/Case: 12 units

Pack Size: 12 x 298g

Packaging Type: Cardboard master

Case Weight Net: 3.576kg

Case Weight Gross: 4.296kg

Case Dimensions: L17.19 x W26.99 x H35.19 cm

Ti/Hi: 11 x 9

Case Cube: 0.01633

Storage Condition: Frozen

Storage Temperature: -18°/0°F

Date Code Format: Best Before YYYY-MM-DD

Shelf Life: 540 days

Valeur nutritive

Nutrition Facts

Pour 1 plat (298 g)

Per 1 tray (298 g)

Calories 420

% valeur quotidienne *
% Daily Value *

Lipides / Fat 19 g	25 %
saturés / Saturated 11 g	58 %
+ trans / trans 0,5 g	

Glucides / Carbohydate 44 g

Fibres / Fibre 2 g	7 %
Sucres / Sugars 2 g	2 %

Protéines / Protein 18 g

Cholestérol / Cholesterol 55 mg

Sodium 700 mg	30 %
----------------------	-------------

Potassium 130 mg	3 %
------------------	------------

Calcium / Calcium 355 mg	27 %
--------------------------	-------------

Fer / Iron 2 mg	9 %
-----------------	------------

* 5% ou moins c'est **peu**, 15% ou plus c'est **beaucoup**

* 5% or less is a **little**, 15% or more is a **lot**

Allergens: **WHEAT, MILK**

Ingredients: Cavatappi pasta [Water, Cavatappi pasta (durum wheat semolina), Cheddar cheese (Milk, bacterial culture, salt, calcium chloride, microbial enzyme, powdered cellulose, natural color), Cream (cream, milk), White old cheddar cheese, Canola oil, Modified corn starch, Skim milk powder, Salt, Garlic, Powdered cheese flavor (yeast extract, salt, natural flavors), Xanthan gum, Spices, Spices extract].

Heating Instructions/Notes: Microwave: Pierce the plastic film with a fork. Heat on high for 4 minutes and mix well. Tray will be hot. Remove from the microwave oven with caution. Let stand for 1 minute before enjoying.

Micro-onde: Percer la pellicule de plastique avec la fourchette. Chauffer pendant 4 minutes et remuer. Retirer avec précaution du micro-onde car le plat sera chaud. Laisser reposer 1 minute. Le temps peut varier selon la puissance des appareils.